



List 1: Shipping solvent of tetrahydrofuran can be replaced with: benzene, chloroform, toluene, xylene, dichloromethane or dischloroethane. Shipping solvent of acetone can be replaced with: carbon tetrachloride, o-dichlorobenzene, dimethylformamide, dodecane, dimethyl sulfoxide, dioxane, ethylacetate, FC-113, hexane, pyridine, hexafluoroisopropanol, chloroform, methyl ethyl ketone, quinolone or cyclohexane. Shipping solvent of chloroform can be replaced with: m-cresol in chloroform or up to 10% hexafluoroisopropanol/chloroform. Shipping solvent of dimethylformamide can be replaced with: dimethyl sulfoxide, dioxane, tetrahydrofuran or toluene. Shipping solvent of o-dichlorobenzene can be replaced with: 1-chloronaphthalene or trichlorobenzene.

List 2: Shipping solvent of tetrahydrofuran can be replaced with: acetone, ethanol, quinolone, benzene, o-dichlorobenzene, ethyl acetate, dodecane, FC-113, carbon tetrachloride, dichloromethane, dichloroethane, n-hexane, cyclohexane, xylene, tetrahydrofuran, chloroform, 1,4-dioxane, hexafluoroisopropanol, toluene, 1-chloronaphthalene, N,N-dimethylacetamide, methyl ethyl ketone, trichlorobenzene, m-cresol, dimethylformamide, methylpyrrolidone, o-chlorophenol/chloroform, dimethyl sulfoxide, and pyridine.

List 3: Shipping solvent of water can be replaced with: aqueous buffers, methanol, ethanol, tetrahydrofuran, dimethylformamide, hexafluoroisopropanol, dimethyl sulfoxide, 2-propanol, and acetonitrile.

List 4: Shipping solvent of water can be replaced with: aqueous buffers with up to 20% methanol, ethanol, propanol, acetonitrile, dimethylformamide, dimethyl sulfoxide, formic acid, and acetic acid and up to 50% acetone.